

## What is Staub Cranial?

The Staub Cranial Technique is the first and only system of its kind. We are able to combine math and science to restore your natural tooth length, width, look and position. By using laws of nature we are able to assist you in providing a fast, accurate and natural looking restoration.

Contact us at Crown Dental Laboratory of Lakeland, and let us demonstrate how the Staub Cranial Technique can be used to help with your dental restoration. The Staub Cranial Technique can be used in just about any restoration, full Dentures, long span bridges, or implant cases.

Crown Dental Lab can save you time and money on your dental restorations. Contact us for more information.

[www.crowndental.biz](http://www.crowndental.biz)

**STAUB**  
CRANIAL™



## Natural Restoration

Combining Math and Science to provide your patients with **dentures in just 3 appointments.**

**Crown Dental Laboratory**  
325 Tyler Avenue  
Lakeland, FL 33801

(863)688-8608  
(800) 860-8608

# The Staub Cranial System

The Staub Cranial System is revolutionary in its thinking, methods, and science. Invented in Neu-Ulm Germany, The Staub Cranial system has exploded in Europe, and is just now making its introduction to North America.

This technology **has the ability to eliminate appointments, saving you time and money.** You are able to get your final restoration back faster with precise results to avoid having to make multiple follow up appointments for adjustments.

Staub Cranial has been proven to give a more accurate, more cosmetic, and longer lasting restoration. Ask us about The Staub Cranial System and start to experience how math and science in dentistry makes the difference.



## How it works –

The Staub Cranial System is a simple but effective way to locate natural tooth position. A restoration of this kind can result in an easier way for your patients to enjoy the simple things. Once again, they will enjoy a full range of motion; eat, laugh, sing, whistle and sneeze.

This system analyzes bone structure landmarks in the mouth. These landmarks are symmetrical, and they work as a key to unlock the natural tooth position. We are able to use these measurements, combined with mathematical and scientific formulas, to provide the best possible restoration.


## Get it done right...the first time.

Below is a picture of multiple removable restorations for the same patient. She was using them all. She had one set for eating, one set for pictures, and one for talking. We were able to replace her bag of dentures with one Staub Cranial made denture. By using exact measurements, we are able to get exact results.



**In Just 3 Appointments!!**

**If you can give us these landmarks in your impression**



**We can give you exact tooth placement.....96.4% of the time....**

